

Anti-IL18 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q14116

Item Number: CM0012928



Introduction

Rabbit polyclonal antibody targeting human IL-18 (IFN-gamma-inducing factor), validated for WB, IP, ELISA. Ideal for studying inflammatory cytokine signaling, Th1/NK cell responses, pyroptosis, and ovarian carcinoma research. - Q14116

[Learn More](#)

Parameter	Value
Product Name	IL18 Rabbit pAb
Remarks/Alias	IGIF; IL-18; IL-1g; IL1F4; IL18
Species	Human
GeneID (Human)	3606
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 37-193 of human IL18 (NP_001553.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	Q14116
Protein Weight	22kDa
Shipping	Ice bag