

Anti-IL18 Rabbit Polymab® For Wb, If/Icc, If-P, Ihc-P, Elisa - Q14116

Item Number: CM0004289

Introduction

High-quality rabbit monoclonal antibody against human IL18 (mature IL-18). Validated for WB, IF/ICC, IF-P, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for studying pro-inflammatory cytokine IL-18 in epithelial barrier, Th1/NK cell responses, and pyroptosis.

[Learn More](#)

Parameter	Value
Product Name	IL18 Rabbit PolymAb®
Remarks/Alias	mature IL-18
Species	Human
GeneID (Human)	3606
GeneID	3606
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-193 of human IL18 (NP_001553.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14116
Protein Weight	18 KDa
Shipping	Ice bag