

Anti-IL18Bp Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O95998

Item Number: CM0028431



Introduction

Rabbit polyclonal antibody against human IL18BP (interleukin-18-binding protein) for WB, IF/ICC, and ELISA. Cross-reacts with mouse IL18BP. Useful for studying IL-18 inhibition and early TH1 cytokine response.

[Learn More](#)

Parameter	Value
Product Name	IL18BP Rabbit pAb
Remarks/Alias	FVH, IL18BPα, IL18BP
Species	Human
Gene ID (Human)	10068
Gene ID	10068
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-194 of human IL18BP (NP_001138529.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95998
Protein Weight	21kDa
Shipping	Ice bag