


Anti-IL18R1/IL-18 R Alpha Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13478

Item Number: CM0023710



Introduction

High-quality rabbit polyclonal antibody targeting human IL18R1 (CD218a, IL-18R α), suitable for Western blot, immunohistochemistry (IHC-P), immunofluorescence (IF/ICC), and ELISA. Detects human, mouse, and rat IL18R1. Ideal for investigating IL18-mediated Th1 cell activation and inflammatory signaling.

[Learn More](#)

Parameter	Value
Product Name	IL18R1/IL-18 R alpha Rabbit pAb
Remarks/Alias	CD218a; IL18RA; IL1RRP; CDw218a; IL-1Rrp; IL-18Ralpha; IL18Ralpha2; IL-18R-alpha; IL18R1/IL-18 R alpha
Species	Human
GeneID (Human)	8809
GeneID	8809
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13478
Protein Weight	62kDa
Shipping	Ice bag