

Anti-IL1R1 Rabbit Polyclonal Antibody For Wb And Elisa - P14778

Item Number: CM0026442



Introduction

Rabbit polyclonal antibody against human IL1R1, validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for research on interleukin-1 signaling and inflammatory response.

[Learn More](#)

Parameter	Value
Product Name	IL1R1 Rabbit pAb
Remarks/Alias	P80; IL1R; IL1RA; CD121A; D2S1473; IL-1RT1; IL-1R-alpha; IL1R1
Species	Human
Gene ID (Human)	3554
Gene ID	3554
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 380-569 of human IL1R1 (NP_000868.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14778
Protein Weight	65kDa
Shipping	Ice bag