


Anti-IL1R2 Rabbit Polyclonal Antibody For Wb, Elisa - P27930

Item Number: CM0023493



Introduction

Rabbit polyclonal antibody against human IL1R2, validated for WB and ELISA. Recognizes both human and mouse IL1R2, the decoy receptor for IL1 beta involved in inflammation and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	IL1R2 Rabbit pAb
Remarks/Alias	IL1RB; CD121b; IL1R2c; CDw121b; IL-1R-2; IL-1RT2; IL-1RT-2; IL1R2
Species	Human
Gene ID (Human)	7850
Gene ID	7850
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-320 of human IL1R2 (NP_775465.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P27930
Protein Weight	45kDa
Shipping	Ice bag