

Anti-IL1Ra Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - P18510

Item Number: CM0022797



Introduction

Rabbit polyclonal antibody against IL1RA (Interleukin-1 receptor antagonist) for Western blot, IF/ICC, IF-P, and ELISA applications. Detects human, mouse, and rat IL1RA. Useful for inflammation and cytokine signaling studies.

[Learn More](#)

Parameter	Value
Product Name	IL1RA Rabbit pAb
Remarks/Alias	DIRA; IRAP; IL1F3; IL1RA; MVCD4; IL-1RN; IL-1ra; IL-1ra3; ICIL-1RA
Species	Human
GeneID (Human)	3557
GeneID	3557
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-177 of human IL1RA (NP_776214.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18510
Protein Weight	20kDa
Shipping	Ice bag