


Anti-IL1RL1 Rabbit Polyclonal Antibody For Wb, Elisa - Q01638

Item Number: CM0022913



Introduction

IL1RL1 Rabbit pAb validated for WB and ELISA; detects human, mouse, rat IL1RL1 (ST2), a receptor for IL-33 involved in innate and adaptive immunity. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	IL1RL1 Rabbit pAb
Remarks/Alias	T1; ST2; DER4; ST2L; ST2V; FIT-1; IL33R; IL1RL1
Species	Human
GeneID (Human)	9173
GeneID	9173
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-328 of human IL1RL1 (NP_057316.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01638
Protein Weight	63kDa
Shipping	Ice bag