


Anti-IL1B Rabbit Monoclonal Antibody For Wb, Elisa - P10749

Item Number: CM0004447



Introduction

IL1 β Rabbit Monoclonal Antibody validated for WB and ELISA. Detects mouse IL1 β (P10749), a pro-inflammatory cytokine involved in fever, T-cell activation, and pyroptosis. Suitable for inflammation and immune research.

[Learn More](#)

Parameter	Value
Product Name	IL1 β Rabbit mAb
Remarks/Alias	IL-1b; IL-1beta; IL1 β
Species	Mouse
GeneID (Human)	3553
GeneID	16176
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 100-269 of mouse IL1 β (P10749).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P10749
Protein Weight	31kDa
Shipping	Ice bag