

Anti-IL1B Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P01584

Item Number: CM0004905



Introduction

High-specificity IL1 β Rabbit mAb validated for WB, IF-P, ELISA. Reacts with human. Targets a pro-inflammatory cytokine key in immune response and disease research. Ideal for IVD development.

[Learn More](#)

Parameter	Description
Product Name	IL1 β Rabbit mAb
Remarks/Alias	IL-1; IL1F2; IL1beta; IL1-BETA; IL1 β
Species	Human
Gene ID (Human)	3553
Gene ID	3553
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-269 of IL1 β (NP_000567.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human
SWISS	P01584
Protein Weight	31 kDa
Shipping	Ice bag