

## Anti-IL1B Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P01584

Item Number: CM0005315



### Introduction

Rabbit monoclonal antibody targeting human IL-1 $\beta$  (interleukin-1 beta), a key pro-inflammatory cytokine. Validated for WB, IF/ICC, and ELISA. Ideal for studying inflammation, pyroptosis, and immune responses.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | IL1 $\beta$ Rabbit mAb                                                                                                   |
| Remarks/Alias    | IL-1; IL1F2; IL1beta; IL1-BETA                                                                                           |
| Species          | Human                                                                                                                    |
| Gene ID (Human)  | 3553                                                                                                                     |
| Gene ID          | 3553                                                                                                                     |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 117-269 of human IL1 $\beta$ (NP_000567.1) |
| Source           | Rabbit                                                                                                                   |
| Category         | Monoclonal Antibodies                                                                                                    |
| Application      | WB, IF/ICC, ELISA                                                                                                        |
| Cross Reactivity | Human                                                                                                                    |
| SWISS            | P01584                                                                                                                   |
| Protein Weight   | 31kDa                                                                                                                    |
| Shipping         | Ice bag                                                                                                                  |