


Anti-IL1B Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P10749

Item Number: CM0005857



Introduction

High-affinity rabbit monoclonal antibody targeting mouse IL1 β , validated for western blot, immunofluorescence, and ELISA. Ideal for studying inflammation, cytokine signaling, and pyroptosis in activated macrophages and related models.

[Learn More](#)

Parameter	Value
Product Name	IL1 β Rabbit mAb
Remarks/Alias	IL-1b; IL-1beta; IL1 β
Species	Mouse
Gene ID (Human)	3553
Gene ID	16176
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 118-269 of mouse IL1 β (NP_032387.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	P10749
Protein Weight	31 kDa
Shipping	Ice bag