


## Anti-IL1B Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P10749

Item Number: CM0008229



### Introduction

Recombinant IL1 $\beta$  rabbit monoclonal antibody validated for WB, IF/ICC, IHC-P, and ELISA. Reacts with human, mouse, and rat IL1 $\beta$ . Ideal for inflammation and cytokine research.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | IL1 $\beta$ Rabbit mAb                                                                                                   |
| Remarks/Alias    | IL-1b; IL-1beta; IL1 $\beta$                                                                                             |
| Species          | Mouse                                                                                                                    |
| GeneID (Human)   | 3553                                                                                                                     |
| GeneID           | 16176                                                                                                                    |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 118-269 of mouse IL1 $\beta$ (NP_032387.1) |
| Source           | Rabbit                                                                                                                   |
| Category         | Monoclonal Antibodies                                                                                                    |
| Application      | WB, IF/ICC, IHC-P, ELISA                                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                                        |
| SWISS            | P10749                                                                                                                   |
| Protein Weight   | 31kDa                                                                                                                    |
| Shipping         | Ice bag                                                                                                                  |