


Anti-IL2 Rabbit Monoclonal Antibody For Wb, Elisa - P04351

Item Number: CM0003221



Introduction

Recombinant rabbit monoclonal antibody against mouse IL2, validated for Western blot and ELISA. Suitable for immune response research, targeting the T-cell growth factor at 19 kDa.

[Learn More](#)

Parameter	Value
Product Name	IL2 Rabbit mAb
Remarks/Alias	IL-2; IL2
Species	Mouse
Gene ID (Human)	3558
Gene ID	16183
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-169 of mouse IL2 (NP_032392.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P04351
Protein Weight	19kDa
Shipping	Ice bag