

Anti-IL2 Rabbit Monoclonal Antibody For Wb, Fc, Elisa - P60568

Item Number: CM0002792



Introduction

High-quality rabbit monoclonal antibody against human Interleukin-2 (IL-2). Suitable for Western blot, flow cytometry, and ELISA. Researchers use it to study immune response, T-cell proliferation, and cytokine signaling in immunotherapy development.

[Learn More](#)

Parameter	Value
Product Name	IL2 Rabbit mAb
Remarks/Alias	IL-2; TCGF; Lymphokine; IL2
Species	Human
Gene ID (Human)	3558
Gene ID	3558
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-153 of human IL2 (NP_000577.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	P60568
Protein Weight	18kDa
Shipping	Ice bag