


Anti-IL2 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P04351

Item Number: CM0007182



Introduction

Recombinant rabbit monoclonal antibody targeting mouse interleukin-2 (IL2). Validated for Western blot, immunofluorescence, and ELISA. Cross-reacts with rat. Ideal for studying immune regulation and T-cell function.

[Learn More](#)

Parameter	Value
Product Name	IL2 Rabbit mAb
Remarks/Alias	IL-2
Species	Mouse
Gene ID (Human)	3558
Gene ID (Mouse)	16183
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 21-169 of mouse IL2 (NP_032392.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P04351
Protein Weight	19kDa
Shipping	Ice bag