


Anti-IL2 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P60568

Item Number: CM0010297



Introduction

Recombinant rabbit monoclonal antibody targeting human IL-2 (interleukin-2, TCGF). Suitable for Western blot, immunofluorescence, and ELISA applications. Ideal for studying T-cell growth, immune response, and tolerance.

[Learn More](#)

Parameter	Specification
Product Name	IL2 Rabbit mAb
Remarks/Alias	IL-2; TCGF; lymphokine
Species	Human
Gene ID (Human)	3558
Gene ID	3558
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-153 of human IL-2 (NP_000577.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P60568
Protein Weight	18kDa
Shipping	Ice bag