

Anti-IL2 Rabbit Polyclonal Antibody For Wb, Elisa - P60568

Item Number: CM0013277



Introduction

Rabbit polyclonal antibody targeting human Interleukin-2 (IL-2). Validated for WB and ELISA applications. Ideal for studying T-cell proliferation, immune tolerance, and cytokine signaling.

[Learn More](#)

Parameter	Value
Product Name	IL2 Rabbit pAb
Remarks/Alias	IL-2; TCGF; lymphokine; IL2
Species	Human
GeneID	3558
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-153 of human IL2 (NP_000577.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P60568
Protein Weight	18 kDa
Shipping	Ice bag