


Anti-IL2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P60568

Item Number: CM0023605



Introduction

High-quality rabbit polyclonal antibody against human IL-2 (IL2), validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for immune response and T-cell research.

[Learn More](#)

Parameter	Value
Product Name	IL2 Rabbit pAb
Remarks/Alias	IL-2; TCGF; Lymphokine; IL2
Species	Human
GeneID (Human)	3558
GeneID	3558
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-153 of IL2 (NP_000577.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60568
Protein Weight	18kDa
Shipping	Ice bag