


Anti-IL20Rb Polyclonal Antibody For Wb, Elisa - Q6Uxl0

Item Number: CM0024819



Introduction

CamelBio offers this rabbit anti-IL20RB polyclonal antibody, validated in WB and ELISA, with cross-reactivity to human, mouse, and rat, for cytokine receptor and inflammatory research.

[Learn More](#)

Parameter	Value
Product Name	IL20RB Rabbit pAb
Remarks/Alias	DIRS1; FNDC6; IL-20R2; IL20RB
Species	Human
Gene ID (Human)	53833
Immunogen	Recombinant fusion protein corresponding to amino acids 48-177 of human IL20RB (NP_653318.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6UXL0
Protein Weight	35kDa
Shipping	Ice bag