

Anti-IL21 Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q9Hbe4

Item Number: CM0030122



Introduction

IL21 Rabbit pAb is a polyclonal antibody validated for WB, IF/ICC, IHC-P, ELISA. It detects human, mouse, and rat IL-21, a key cytokine in immune regulation. Suitable for research on Tfh cells, germinal centers, and antibody responses.

[Learn More](#)

Parameter	Value
Product Name	IL21 Rabbit pAb
Remarks/Alias	Za11; IL-21; CVID11; IL21
Species	Human
Gene ID	59067
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-162 of human IL21 (NP_068575.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HBE4
Protein Weight	19 kDa
Shipping	Ice bag