


Anti-IL22 Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Gzx6

Item Number: CM0004307



Introduction

Rabbit monoclonal antibody against human IL22 (Interleukin-22), validated for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for inflammation and epithelial barrier research.

[Learn More](#)

Parameter	Value
Product Name	IL22 Rabbit mAb
Remarks/Alias	TIFa; IL-21; IL-22; ILTIF; IL-TIF; IL-D110; zcyto18; TIFIL-23; IL22
Species	Human
Gene ID (Human)	50616
Gene ID	50616
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 34-179 of human IL22 (NP_065386.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZX6
Protein Weight	20kDa
Shipping	Ice bag