

Anti-IL22 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Jjy9

Item Number: CM0018793



Introduction

Rabbit polyclonal antibody targeting mouse IL22, validated for WB and ELISA. Cross-reacts with human and mouse. Recognizes a cytokine essential for epithelial barrier integrity and tissue regeneration. - Q9JJY9

[Learn More](#)

Parameter	Value
Product Name	IL22 Rabbit pAb
Remarks/Alias	IL-22; If2b1; Iltif; IL-22a; ILTIFa; IL22
Species	Mouse
Gene ID	50929
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 34-179 of mouse IL22. (NP_058667.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9JJY9
Protein Weight	20kDa
Shipping	Ice bag