


Anti-IL22 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Gzx6

Item Number: CM0029871



Introduction

Rabbit polyclonal antibody against human IL22, suitable for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying inflammation, epithelial barrier function, and regenerative medicine. SWISS: Q9GZX6.

[Learn More](#)

Parameter	Value
Product Name	IL22 Rabbit pAb
Remarks/Alias	TIFa; IL-21; IL-22; ILTIF; IL-TIF; IL-D110; zcyto18; TIFIL-23; IL22
Species	Human
Gene ID (Human)	50616
Gene ID	50616
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 34-179 of human IL22 (NP_065386.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZX6
Protein Weight	20kDa
Shipping	Ice bag