

Anti-IL22Ra1 Polyclonal Antibody For Wb, Elisa - Q8N6P7

Item Number: CM0016181



Introduction

Rabbit polyclonal antibody against human IL22RA1 (IL22R, CRF2-9), suitable for WB and ELISA. Cross-reacts with mouse and rat. Predicted molecular weight 63 kDa. A valuable tool for IL-22 signaling research.

[Learn More](#)

Parameter	Value
Product Name	IL22RA1 Rabbit pAb
Remarks/Alias	IL22R; CRF2-9; IL22R1; IL22RA1
Species	Human
Gene ID (Human)	58985
Gene ID	58985
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 16-230 of human IL22RA1 (NP_067081.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N6P7
Protein Weight	63kDa
Shipping	Ice bag