


Anti-IL23R Monoclonal Antibody For Fc - Q5Vwk5

Item Number: CM0011771



Introduction

Rabbit monoclonal antibody against human interleukin-23 receptor (IL23R), conjugated to . Designed for flow cytometry (FC). Recognizes the receptor involved in T-cell and NK cell activation, with implications in autoimmunity and inflammation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IL23R mAb
Remarks/Alias	PSORS7
Species	Human
GeneID (Human)	149233
GeneID	149233
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q5VWK5
Protein Weight	41kDa/43kDa/72kDa
Shipping	Ice bag