


Anti-IL23R Rabbit Monoclonal Antibody For Fc, Elisa - Q5Vwk5

Item Number: CM0009010



Introduction

Anti-IL23R Rabbit Monoclonal Antibody for flow cytometry (FC) and ELISA. Recognizes human IL23R (SwissProt Q5VWK5). A key tool for studying IL-23 signaling in innate and adaptive immunity.

[Learn More](#)

Parameter	Value
Product Name	IL23R Rabbit mAb
Remarks/Alias	PSORS7
Species	Human
GeneID (Human)	149233
GeneID	149233
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 124-353 of human IL23R (NP_653302.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC, ELISA
Cross Reactivity	Human
SWISS	Q5VWK5
Protein Weight	41kDa/43kDa/72kDa
Shipping	Ice bag