


Anti-IL23R Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q5Vwk5

Item Number: CM0003320



Introduction

High-quality rabbit monoclonal antibody against human IL23R (alias PSORS7). Validated for WB, IHC-P, and ELISA with cross-reactivity to human, mouse, and rat. Ideal for IL23 receptor signaling and immunology research.

[Learn More](#)

Parameter	Value
Product Name	IL23R Rabbit mAb
Remarks/Alias	PSORS7; IL23R
Species	Human
GeneID	149233
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 510-629 of human IL23R (NP_653302.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5VWK5
Protein Weight	72kDa
Shipping	Ice bag