


## Anti-IL23R Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q5Vwk5

Item Number: CM0021148



### Introduction

Rabbit polyclonal antibody targeting human IL23R, suitable for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Immunogen: recombinant fusion protein (aa 510-629). Ideal for IL23/Th17 pathway and autoimmune research.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | IL23R Rabbit pAb                                                                                                    |
| Remarks/Alias    | PSORS7; IL23R                                                                                                       |
| Species          | Human                                                                                                               |
| Gene ID (Human)  | 149233                                                                                                              |
| Gene ID          | 149233                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 510-629 of human IL23R (NP_653302.2). |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, IF/ICC, IHC-P, ELISA                                                                                            |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | Q5VWK5                                                                                                              |
| Protein Weight   | 72kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |