


Anti-IL24 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q13007

Item Number: CM0020188



Introduction

IL24 rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. For IL24-mediated immune response, oncogenesis, and antiviral cytokine signaling research.

[Learn More](#)

Parameter	Value
Product Name	IL24 Rabbit pAb
Remarks/Alias	C49A; FISP; MDA7; MOB5; ST16; IL10B; IL24
Species	Human
GeneID (Human)	11009
GeneID	11009
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 52-140 of human IL24 (NP_001172085.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q13007
Protein Weight	24kDa
Shipping	Ice bag