


Anti-IL25 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q9H293

Item Number: CM0028792



Introduction

Rabbit polyclonal antibody against human IL25 (IL-17E). Validated for WB, IF-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for investigating Th2 immune responses and cytokine signaling pathways.

[Learn More](#)

Parameter	Specification
Product Name	IL25 Rabbit pAb
Remarks/Alias	IL17E; IL25
Species	Human
Gene ID (Human)	64806
Gene ID	64806
Immunogen	Recombinant fusion protein containing amino acids 33-177 of human IL25 (NP_073626.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H293
Protein Weight	20kDa
Shipping	Ice bag