

Anti-IL26 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Nph9

Item Number: CM0005754

Introduction

High-quality rabbit monoclonal antibody specific for human Interleukin-26 (IL26). Validated for Western blotting and ELISA. Suitable for mucosal immunity and inflammatory bowel disease research. Predicted band ~20kDa.

[Learn More](#)

Parameter	Value
Product Name	IL26 Rabbit mAb
Remarks/Alias	AK155; IL-26; IL26
Species	Human
GeneID (Human)	55801
GeneID	55801
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-154 of human IL26 (NP_060872.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NPH9
Protein Weight	20kDa
Shipping	Ice bag