

Anti-IL26 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nph9

Item Number: CM0021753

Introduction

High-quality rabbit polyclonal antibody targeting human IL26, validated for Western blotting (WB) and ELISA. Cross-reacts with mouse. Suitable for mucosal immunity, inflammatory bowel disease, and cytokine signaling studies.

[Learn More](#)

Parameter	Value
Product Name	IL26 Rabbit pAb
Remarks/Alias	AK155; IL-26; IL26
Species	Human
GeneID (Human)	55801
GeneID	55801
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-100 of human IL26 (NP_060872.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9NPH9
Protein Weight	20kDa
Shipping	Ice bag