


Anti-IL2Rb (Cd122/IL15Rb) Polyclonal Antibody For Wb, Elisa - P14784

Item Number: CM0028232



Introduction

High-quality rabbit polyclonal antibody targeting human IL2RB (CD122/IL15RB), validated for Western blot (WB) and ELISA. Shows cross-reactivity with mouse and rat. Suitable for IL-2 and IL-15 receptor signaling studies.

[Learn More](#)

Parameter	Value
Product Name	IL2RB Rabbit pAb
Remarks/Alias	CD122; IMD63; IL15RB; P70-75; IL2RB
Species	Human
Gene ID (Human)	3560
Gene ID	3560
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 316-431 of human IL2RB (NP_000869.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14784
Protein Weight	61kDa
Shipping	Ice bag