

Anti-IL3 Rabbit Polyclonal Antibody For Wb, Elisa - P08700

Item Number: CM0030912



Introduction

Rabbit polyclonal antibody targeting IL-3 (interleukin-3, UniProt P08700), validated for WB and ELISA. Cross-reacts with human and mouse IL-3, a key hematopoietic growth factor. Recombinant immunogen (aa20-152). Ideal for cytokine and immunology research.

[Learn More](#)

Parameter	Value
Product Name	IL3 Rabbit pAb
Remarks/Alias	IL-3; MCGF; MULTI-CSF; IL3
Species	Human
GeneID (Human)	3562
GeneID	3562
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-152 of human IL3 (NP_000579.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P08700
Protein Weight	17kDa
Shipping	Ice bag