

Anti-IL31 Rabbit Polyclonal Antibody For Wb, Elisa - Q6Ebc2

Item Number: CM0015777

Introduction

Rabbit polyclonal antibody against human IL-31, validated for Western blot and ELISA. Recognizes human, mouse, and rat IL-31. Designed for research on cytokine signaling and skin immunity.

[Learn More](#)

Parameter	Value
Product Name	IL31 Rabbit pAb
Remarks/Alias	IL-31; IL31
Species	Human
Gene ID (Human)	386653
Gene ID	386653
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-164 of human IL31 (NP_001014358.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6EBC2
Protein Weight	18kDa
Shipping	Ice bag