

Anti-IL31Ra Polyclonal Antibody For Wb, Elisa - Q8Ni17

Item Number: CM0029612



Introduction

Rabbit polyclonal antibody against human Interleukin-31 receptor subunit alpha (IL31RA). Suitable for Western blot and ELISA, cross-reactive with mouse. Explore IL31-mediated itch and skin immunity.

[Learn More](#)

Parameter	Value
Product Name	IL31RA Rabbit pAb
Remarks/Alias	CRL; GPL; CRL3; GLMR; GLM-R; PLCA2; hGLM-R; IL-31RA; PRO21384; zcytoR17; IL31RA
Species	Human
Gene ID	133396 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-321 of human IL31RA (NP_001284501.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8NI17
Protein Weight	83kDa
Shipping	Ice bag