

Anti-IL33 Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - O95760

Item Number: CM0006954



Introduction

Recombinant monoclonal rabbit antibody against human IL-33 (IL-1F11, NF-HEV). Suitable for WB, IF/ICC, IF-P, IHC-P, and ELISA applications. Designed for researchers studying Th2 immune responses, allergic inflammation, and nuclear transcriptional regulation. Ships on ice bag.

[Learn More](#)

Parameter	Value
Product Name	IL33 Rabbit mAb
Remarks/Alias	DVS27; IL1F11; NF-HEV; NFEHEV; C9orf26; IL33
Species	Human
Gene ID (Human)	90865
Gene ID	90865
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 109-270 of human IL33 (NP_001300974.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	O95760
Protein Weight	31kDa
Shipping	Ice bag