


Anti-IL33 Rabbit Polyclonal Antibody For Wb, Elisa - 095760

Item Number: CM0031073



Introduction

CamelBio anti-IL33 rabbit polyclonal antibody for WB and ELISA, detecting human, mouse, and rat IL33. Suitable for studying Th2 differentiation, alarmin signaling, and inflammation.

[Learn More](#)

Parameter	Value
Product Name	Anti-IL33 Rabbit Polyclonal Antibody
Remarks/Alias	DVS27; IL1F11; NF-HEV; NFEHEV; C9orf26; IL33
Species	Human
GeneID (human)	90865
GeneID	90865
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 110-190 of human IL33 (NP_254274.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95760
Protein Weight	31kDa
Shipping	Ice bag