


Anti-IL33 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Bvz5

Item Number: CM0019299



Introduction

Rabbit polyclonal antibody against mouse IL-33, validated for WB and ELISA. Detects mouse and rat IL-33. Suitable for studying Th2 response, alarmin function, and nuclear transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	IL33 Rabbit pAb
Remarks/Alias	IL-33; IL1f11; NF-HEV; 9230117N10Rik; IL-33
Species	Mouse
GeneID (Human)	90865
GeneID (Mouse)	77125
Immunogen	Recombinant fusion protein containing amino acids 109-266 of mouse IL-33 (NP_598536.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8BVZ5
Protein Weight	30kDa
Shipping	Ice bag