

Anti-IL36A Polyclonal Antibody For Wb, Elisa - Q9Uha7

Item Number: CM0030883



Introduction

Rabbit polyclonal antibody against human IL36A (FIL1 epsilon, IL-1F6) for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying IL-36 signaling in inflammatory skin diseases.

[Learn More](#)

Parameter	Value
Product Name	IL36A Rabbit pAb
Remarks/Alias	FIL1; FIL1E; IL1F6; IL-1F6; IL1(EPSILON); FIL1(EPSILON); IL36A
Species	Human
GeneID (Human)	27179
GeneID	27179
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-158 of human IL36A (NP_055255.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UHA7
Protein Weight	18kDa
Shipping	Ice bag