


Anti-IL36G Polyclonal Antibody For Wb, Elisa - Q9Nzh8

Item Number: CM0025748



Introduction

Rabbit polyclonal antibody targeting human IL36G, a cytokine in the IL-1 family. Validated for WB and ELISA; cross-reacts with mouse and rat. Suitable for studying inflammatory responses, psoriasis, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	IL36G Rabbit pAb
Remarks/Alias	IL1E, IL1F9, IL1H1, IL-1F9, IL-1H1, IL1RP2, IL-1RP2, IL36G
Species	Human
GeneID (Human)	56300
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-169 of human IL36G (NP_062564.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZH8
Protein Weight	19kDa
Shipping	Ice bag