

Anti-IL36R Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Hb29

Item Number: CM0014708



Introduction

Rabbit polyclonal antibody against human IL36R (IL-1 receptor-like 2), validated for WB, IF/ICC, ELISA. Reacts with human and mouse. Suitable for studying IL-36 signaling in inflammatory responses.

[Learn More](#)

Parameter	Value
Product Name	IL36R Rabbit pAb
Remarks/Alias	IL-36R; IL1RRP2; IL-1Rrp2; IL1R-rp2; IL36R
Species	Human
Gene ID (Human)	8808
Gene ID	8808
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-116 of human IL36R (NP_003845.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9HB29
Protein Weight	65kDa
Shipping	Ice bag