


Anti-IL3Ra Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P26951

Item Number: CM0025031



Introduction

Rabbit polyclonal antibody targeting human IL3RA (CD123), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for hematopoietic cell and leukemia studies.

[Learn More](#)

Parameter	Value
Product Name	IL3RA Rabbit pAb
Remarks/Alias	IL3R, CD123, IL3RX, IL3RY, IL3RAY, hIL-3Ra, IL-3R-alpha, IL3RA
Species	Human
Gene ID (Human)	3563
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human IL3RA (NP_002174.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26951
Protein Weight	43kDa
Shipping	Ice bag