


Anti-IL4 Rabbit Monoclonal Antibody For Wb And Elisa - P05112

Item Number: CM0005739



Introduction

Recombinant rabbit monoclonal antibody targeting human interleukin-4 (IL4), validated for Western blot and ELISA applications. Ideal for investigating B-cell activation, IgE class switching, and inflammatory cytokine signaling.

[Learn More](#)

Parameter	Value
Product Name	IL4 Rabbit mAb
Remarks/Alias	BSF1; IL-4; BCGF1; BSF-1; BCGF-1; IL4
Species	Human
Gene ID (Human)	3565
Gene ID	3565
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-153 of human IL4 (NP_000580.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P05112
Protein Weight	17kDa
Shipping	Ice bag