

Anti-IL4 Rabbit Monoclonal Antibody For Wb, Elisa - P07750

Item Number: CM0006910

Introduction

Rabbit monoclonal antibody targeting mouse IL-4, validated for WB and ELISA. Ideal for immune regulation, inflammation, and B-cell response research. - P07750.

[Learn More](#)

Parameter	Value
Product Name	IL4 Rabbit mAb
Remarks/Alias	IL-4; BSF-1; IL4
Species	Human
GeneID (Human)	3565
GeneID	16189
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-140 of mouse IL4 (NP_067258.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P07750
Protein Weight	16kDa
Shipping	Ice bag