

## Anti-IL4I1 Polyclonal Antibody For Wb, Elisa - Q96Rq9

Item Number: CM0025950

### Introduction

Rabbit polyclonal antibody against human IL4I1, validated for WB and ELISA. Recognizes human, mouse, and rat IL4I1. Useful for immunology and oncology research involving L-amino-acid oxidase and AHR pathway.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | IL4I1 Rabbit pAb                                                                                                  |
| Remarks/Alias    | LAO; FIG1; LAAO; hIL4I1; IL4I1                                                                                    |
| Species          | Human                                                                                                             |
| GeneID (human)   | 259307                                                                                                            |
| GeneID           | 259307                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 22-220 of human IL4I1 (NP_690863.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | Q96RQ9                                                                                                            |
| Protein Weight   | 63kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |