

Anti-IL4R Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P24394

Item Number: CM0017420



Introduction

Rabbit polyclonal antibody targeting human IL4R (CD124). Validated for WB, IHC-P, ELISA with cross-reactivity to human, mouse, and rat. Suitable for allergy and Th2 differentiation research.

[Learn More](#)

Parameter	Value
Product Name	IL4R Rabbit pAb
Remarks/Alias	CD124; IL4RA; IL-4RA; IL4R
Species	Human
GeneID (Human)	3566
GeneID	3566
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-232 of human IL4R (NP_000409.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24394
Protein Weight	90kDa
Shipping	Ice bag