

Anti-IL5Ra Rabbit Polyclonal Antibody For Wb, Elisa - Q01344

Item Number: CM0030078



Introduction

IL5RA Rabbit pAb validated for WB and ELISA; reacts with human, mouse, rat. Targets interleukin-5 receptor alpha, essential for eosinophil biology. Ideal for IVD assay development and research.

[Learn More](#)

Parameter	Value
Product Name	IL5RA Rabbit pAb
Remarks/Alias	IL5R; CD125; CDw125; HSIL5R3; IL5RA
Species	Human
Gene ID	3568
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-335 of human IL5RA (NP_783855.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01344
Protein Weight	48kDa
Shipping	Ice bag