


Anti-IL6R Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P08887

Item Number: CM0024101



Introduction

Rabbit polyclonal antibody against human IL6R, validated for WB, IHC-P, and ELISA; cross-reacts with mouse and rat. Ideal for IL6 signaling research in inflammation, immunity, and chronic disease.

[Learn More](#)

Parameter	Value
Product Name	IL6R Rabbit pAb
Remarks/Alias	IL6Q; gp80; CD126; HIES5; IL-6R; IL6RA; IL6RQ; IL-1Ra; IL-6RA; IL6QTL; IL-6R-1; IL6R
Species	Human
GeneID	3570
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-320 of human IL6R (NP_000556.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08887
Protein Weight	52kDa
Shipping	Ice bag